
PSI-CAMEA - MAGNET OPTICS

MAGNET OPTICS

Jakob Lass



MARCH 20, 2018
MAGNET OPTICS
NIELS BOHR INSTITUTE
PAUL SCHERRER INSTITUTE

Contents

Table of Contents	A
1 Plan	1
1.1 Notes from meeting 13/03/18	1
1.2 Magnet limitations	1
1.3 Original setup	2
1.4 Comparison	3
1.5 Future Monochromator	5
1.6 Next steps	8
1.7 Folder overview	10

1 Plan

1.1 Notes from meeting 13/03/18

Focusig, in principle for elliptical → simulated by straight for MA15 (test case for MA9)

- 2.3 degrees opening angle 4AA
- $m=6$
- sample size: 5x5mm
- 30cm at the magnet. Compare with only sample position (should be 5x5mm)
- Inner dimation is 7mm (focus for 5 mm on sample position)
- limits for outer diameter
- Reference sample to look at resolution
- ? Maybe Union
- Build up a 'complete' back-end

1.2 Magnet limitations




Figure 1: Sketch of optic dimensions due to magnet limitations.

- MagnetOpticsOpening = 0.030,
- MagnetOpticsEnd = 0.007,
- MagnetOpticsLenght=0.268,
- MagnetOpticsM = 6,
- MagnetOpticsWidth = 0.2,

1.3 Original setup

Setup consists of:

```

COMPONENT MonoFlat = Monochromator_curved(
    zwidth = 0.05, yheight = 0.025, gap = 0.0005, NH = 3,
    NV = 5, DM=3.355, mosaich = 37, mosaicv = 37, r0 = 1,
    RV = RV_mono, RH = RH_mono)
AT (0, 0, 1.6) RELATIVE SlitGuideEnd
ROTATED (0, thetaA1, 0) RELATIVE SlitGuideEnd

COMPONENT A2 = Arm()
AT (0, 0, 1.6) RELATIVE SlitGuideEnd
ROTATED (0, thetaA2, 0) RELATIVE SlitGuideEnd

COMPONENT SlitMagnetStart = Slit(
    xwidth = MagnetSlitWidth, yheight = MagnetSlitHeight)
AT (0, 0, 1.3) RELATIVE A2

COMPONENT PsdMagnet = PSD_monitor(
    nx = 100, ny = 100, filename = "PSD_Magnet_opening.dat", xwidth = 0.2,
    yheight = 0.2)
AT (0, 0, 1.3) RELATIVE A2

COMPONENT DivMon_Magnet = Divergence_monitor(
    nh = 100, nv = 100, filename = "Div_Magnet_opening.dat",
    xwidth = 0.2, yheight = 0.2, maxdiv_h = 6, maxdiv_v = 6.5)
AT (0, 0, 0 ) RELATIVE PREVIOUS

COMPONENT PsdSampleBeamProfile = PSD_monitor(
    nx = 100, ny = 100, filename = "PsdSampleBeamProfile.dat", xwidth = 0.2,
    yheight = 0.2)
AT (0, 0, 1.6) RELATIVE A2

COMPONENT DivMonSampleBeamProfile = Divergence_monitor(
    nh = 100, nv = 100, filename = "DivMonSampleBeamProfile.dat",
    xwidth = 0.2, yheight = 0.2, maxdiv_h = 6, maxdiv_v = 6.5)
AT (0, 0, 0 ) RELATIVE PREVIOUS

COMPONENT PsdSampleSimulator = PSD_monitor(
    nx = 100, ny = 100, filename = "PsdSampleSimulator.dat", xwidth = 0.005,
    yheight = 0.005)
AT (0, 0, 0) RELATIVE PREVIOUS

COMPONENT DivMonSampleSimulator = Divergence_monitor(
    nh = 100, nv = 100, filename = "DivMonSampleSimulator.dat",

```

```
xwidth = 0.005, yheight = 0.005, maxdiv_h = 6, maxdiv_v = 6.5)
AT (0, 0, 0 ) RELATIVE PREVIOUS
```




Figure 2: Flat monochromator with slits and without optics.




Figure 3: Flat monochromator with slits and with optics.

1.4 Comparison

	Flat Mono WO	VCurved Mono WO	Flat Mono W	VCurved Mono W
Intensities	177282	432272	204688	591671
Flat Mono WO	1	2.44	1.15	3.34
VCurved Mono WO	–	1	0.47	1.37
Flat Mono W	–	–	1	2.89
VCurved Mono W	–	–	–	1

Table 1: Comparison between vertically focused/unfocused and with/without optics




Figure 4: No optics flat.




Figure 5: No optics curved vertically.




Figure 6: Flat optics flat.




Figure 7: Flat optics curved vertically.


1.5 Future Monochromator

The future monochromator might change the efficiency of adding the optics. The monochromator is given as

```
COMPONENT Mono = Monochromator_curved(
  zwidth = 0.01, yheight = 0.01, gap = 0.0005, NH =24,
  NV = 13, DM=3.355,mosaich = 37, mosaicv = 37, r0 = 1,
  RV = RV_mono, RH = RH_mono)
AT (0, 0, 1.6) RELATIVE SlitGuideEnd
ROTATED (0, thetaA1,0) RELATIVE SlitGuideEnd
```

	Flat Mono WO	Curved Mono WO	Flat Mono W	VCurved Mono W
Intensities	160212	598328	197003	833920
Flat Mono WO	1	3.73	1.23	5.21
VCurved Mono WO	–	1	0.33	1.39
Flat Mono W	–	–	1	4.23
VCurved Mono W	–	–	–	1

Table 2: Comparison between double focused/unfocused and with/without optics with the future monochromator.






Figure 8: No optics flat with new monochromator.




Figure 9: No optics double focusing with new monochromator.




Figure 10: Flat optics flat with the future monochromator.




Figure 11: Flat optics double focusing with the future monochromator.

1.6 Next steps

Possible next steps

- Elliptical optics
- Optimization of optic opening
- Full simulation



Figure 12: Energy scan, flat monochromator, without optics, VanSample with radius=0.005m and yheight=0.01m




Figure 13: Instrument view of energy scan simulation with two slits.

1.7 Folder overview

1.7.1 Old Monochromator

Flat monochromator without optics:

- OriginalFlat

Vertically curved monochromator without optics:

- OriginalCurved

Flat monochromator with flat optics:

- FlatOpticsFlat

Vertically curved monochromator with flat optics:

- FlatOpticsCurved

1.7.2 New Monochromator

Flat monochromator without optics:

- FutureMonochromator/Flat

Double focusing monochromator without optics:

- FutureMonochromator/DoubleCurved

Flat monochromator with flat optics:

- FutureMonochromator/FlatWithOptics

Doble focusing monochromator with flat optics:

- FutureMonochromator/CurvedWithOptics